# Instruction for use

**Notice:** -Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Intended purpose

Brasseler EZX Blade are medical device for cutting bone , tissue, etc.

## Claims for clinical performance

Cutting bone, being connected with AC-powered, battery-powered, or air-powered devices.

## Contraindication


It is not indicated for long-term use or implantation.

## Complications

Complications related to Brasseler EZX Blade are very rare, but when they occur, they can be serious.

- Patient burn or thermal necrosis
- Paralysis by wrong operation
- Tear or puncture dura or other central nervous system
- Frostbite, hematoma in target site etc. by infection
- Injury by wrong operation or cracked chip of blade

## Explanation of symbols

 Catalogue number	 Batch code	 Medical device
 Caution	 Use by Date	 Keep away from sunlight
 Unique device identifier	 Do not re-use	 Keep dry
 Sterilized using irradiation	 Do not use if package is damaged and consult instructions for use	 Federal law restricts this device to sale by or on the order of a physician.
 Distributor	 Consult instructions for use	 Country of manufacture (Date of manufacture)
 Manufacturer	 Single sterile barrier system with protective packaging inside	
 Do not re-sterilize	 Authorized representative in the European Community/European Union	

## Proper use

- Before use
  - Read Instruction for Use of the product to understand the method of use and precautions.
  - Identify the status of Sterilization.
  - Remove any screws and clear away any soft tissue to ensure no interference with the blade.
- Size the Liner and the Blade
  - Remove the existed acetabular liner.
  - Measure the inner diameter of the cup by using the Liner sizing tool included in the tray.
  - Select the 1st Liner of the same size as the inner diameter of the cup.
  - Place the 1st Liner in the cup.
  - Measure the outer diameter of the cup by using the Blade sizing tool included in the tray.
  - Select the Blade which is 2~3mm larger than the outer diameter of the cup.

## Cooling

Sufficient cooling should be provided when using Brasseler EZX Blade. In case of insufficient cooling, the Brasseler EZX Blade teeth will get clogged with bone chips, leading to increased heat generation. At the worst, this may cause irreversible damage to the bone (thermal necrosis) and reduce service life of the saw blades considerably.

## Disposal

It shall be disposed or managed as an infectious medical waste. It shall be inserted into container(hospital box) with inert double plastic bag(strong, leak proof, puncture proof) after use.

## Caution

- Brasseler EZX Blade are to sale by or on the order of a physician.
- See instruction for use. Important information.
- See detailed handpiece user instruction manual.
- Check the package before use. If package is opened or damaged, do not use.
- Single patient use only.
- Make sure blade is locked and secure before running handpiece and protect the vibrations by handpiece.
- Monitor blade and distal handpiece temperature.
- Irrigate surgical site to reduce heat.
- Do not use excessive force.
- Blunt and damaged Brasseler EZX Blade must no longer be used.
- Do not use twice. This is single use product.
- Storage conditions in dust-free, cool and away from direct sunlight .
- The user of the product is responsible for proper use the equipment during surgery.

- This is single use device not intended for re-sterilization and reuse. The possible outcomes from the reuse of this could be, but not limited to; insufficient or improper cleaning/ registration resulting in cross-infection of patients or degradation of physical characteristics of the device which could pose a risk to the patient. Loss of structural integrity from repeated use is also possible. Detailed examples of possible outcomes from the reuse would be as below:
- Reprocessing a single-use device may alter its characteristics so that it no longer complies with the original manufacturer's specifications and performance may be compromised.
- This device has not undergone extensive testing, validation and documentation to ensure the devices are safe to reuse. Infection is one of the greatest patient safety concerns associated with reuse. The risk of cross infection may increase due to the inability of the reprocessing system to completely remove viable micro-organisms. This may be due to design e.g. narrow lumens. Viable micro-organisms may be incompletely removed and be transferred to the next patient.
- Some materials can absorb certain chemicals which can gradually leach from the material over time. For e.g. disinfectants may be absorbed by plastics and leach out during use, resulting in chemical burns or a risk of sensitization of the patient or user.
- Chemicals may cause corrosion or changes to the materials of the device. For e.g., plastics may soft, crack or become brittle.
- The material may experience stress during reuse and may fail, stretch or break.
- Inadequately cleaned device can carry bacterial endotoxins which remain after bacteria are killed.
- Don't re-start the operation when the blade is staying(remaining) in the bone or tissue, etc. Blade can be broken by unexpected force.

**Brasseler U.S.A Medical, LLC**  
 One Brasseler Boulevard, Savannah, Georgia 31419 U.S.A

**IMEDICOM Co., Ltd.**  
 #612, #709, #804, #815, 172, LS-ro, Gunpo-si, Gyeonggi-do, 15807, Korea  
 Tel : +82 31 479 1156 / Fax : +82 31 479 1157

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## Revendications relatives aux performances cliniques.

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