



BRASSELER
USA[®] SURGICAL

CASE STUDY

Acetabular Revision & Cup Removal in Osteoporotic Bone Using the EZX System



Matthew Deren, MD

Orthopedic Surgery,
Cleveland, Ohio


Patient History & Preoperative Condition:

- Patient: Female with direct anterior total hip arthroplasty performed for femoral neck fracture.
- Acetabular component malpositioned (>30 degrees anteversion).
- Four anterior dislocations following initial surgery.
- Revision surgery planned for acetabular and head components.

Procedure Details:

- Acetabular revision performed due to malpositioning.
- Trunnion impingement on the posterior wall was addressed using a bone hook to clear space for the EZX.
- The EZX system removed the acetabular component in approximately 60 seconds.
- The patient's bone quality was very osteoporotic, but the removal was successful without additional bone loss or fractures.
- The multi-hole acetabular component was 8mm larger than the explanted acetabulum, and a reamer was used to medialize for better coverage.

Pre-Op Image:



Outcomes:

- Device Performance:
 - Provided greater control and confidence.
 - Reduced operative time.
 - Minimal bone loss observed during removal.
- Post-operative Observations:
 - Excellent removal outcome with no additional fractures or complications.
 - The patient had an excellent outcome post-surgery.

Post-Op Image:





BRASSELER

USA® SURGICAL

Bobby Hale

Global Sales Director, Revision Portfolio
(719) 237-4086, BobbyHale@BrasselerUSA.com



Brian Young

Director of Sales, Revision Portfolio West
(214) 405-2962, BrianYoung@BrasselerUSA.com




Steven Cook

Director of Sales, Revision Portfolio East
(919) 538-2206, StevenCook@BrasselerUSA.com



George Blumenschein

Revision Sales Specialist
(973) 487-9642, GeorgeBlumenschein@BrasselerUSA.com



For More Information
Visit us online at **BrasselerSurgical.com**
Call **800.535.6638** to order today!