



**BRASSELER**  
USA® THE SMART CHOICE

## EXTRA-SHARP K-WIRES

Are you concerned with splintering of small bones while using a K-Wire for internal fixation? Try our new extra-sharp K-Wires.

Fabricated from precision ground 316L implant-grade stainless steel, our extra-sharp K-Wires feature a 50° angle point as compared to our standard K-wire with a 90° angle point.

Only available in 4" length, sterile and smooth/plain.

### ► OUR NEW EXTRA-SHARP 50° ANGLE DOUBLE DIAMOND POINT

Diameter		Order Number
mm	in	Sterile (1 Per Pack)
0.7	0.028	KM172-44-28ESS
0.9	0.035	KM172-44-35ESS
1.1	0.045	KM172-44-45ESS
1.6	0.062	KM172-44-62ESS

50° angle point



90° angle point



Call 800-535-6638 to order today!

