



BRASSELER
USA® SURGICAL INSTRUMENTATION

ORTHOPAEDIC PIN & DRILL PACKS


Fabricated from 17-4 stainless steel and are subjected to strenuous inspection programs to assure exceptional results



Stryker® Replacements

Compatible with Stryker® Knee Instrumentation Systems

Order Number	Replacement	Diameter (mm)	Length (mm)	Description	Packaging
BMA-S90008S	7650-2038 7650-1038	3.1	88.9	Headless Pin (Square End & Trocar Point)	4 Pins Per Pack Sterile
BMA-S90010S	7650-2038A 7650-1038A	3.1	88.9	Headless Pin, Fluted (Square End & Trocar Point)	4 Pins Per Pack Sterile
BMA-S90012S	6541-4-003 6541-004-003	3.1	76.2	Headless Pin (Square End & Trocar Point)	4 Pins Per Pack Sterile



Compatible with Stryker® Ortholock™ Navigation System

Order Number	Replacement	Diameter (mm)	Length (mm)	Description	Packaging
ST6007-003-100	6007-003-100	3	100	Navigation Pin (Square End & Trocar Point)	2 Pins Per Pack Sterile
ST6007-003-150	6007-003-150	3	150	Navigation Pin (Square End & Trocar Point)	2 Pins Per Pack Sterile
ST6007-103-110	6007-103-110	3	110	Navigation Pin (Square End & Trocar Point)	2 Pins Per Pack Sterile
ST6007-104-150	6007-104-150	4	150	Navigation Pin (Square End & Trocar Point)	2 Pins Per Pack Sterile

Stryker® and OrthoLock™ are trademarks of Stryker Corporation or its affiliates.

Call 800.535.6638 to order today! Visit us online at BrasselerUSAMedical.com





BRASSELER
USA® SURGICAL INSTRUMENTATION

ORTHOPAEDIC PIN & DRILL PACKS

Fabricated from 17-4 stainless steel and are subjected to strenuous inspection programs to assure exceptional results



DePuy® Synthes® Replacements

Compatible with DePuy® Synthes® Knee Instrumentation Systems

Order Number	Replacement	Diameter (mm)	Length (mm)	Description	Packaging
BMA-S90006S	2490-95-000 2490-95	3.1	76.2	Fixation Pin (Round End & Trocar Point)	4 Pins Per Pack Sterile
BMA-S90007S	86-4192-000 86-4192	3.1	76.2	Fixation Pin (Triangular End & Trocar Point) & Twist Drills (Round End)	4 Pins & 2 Drills Per Pack Sterile
DP86-9117	86-9117 869117	3.1	76.2	Fixation Pin (Triangular End & Trocar Point) & Twist Drills (Round End)	10 Pins Per Pack Non-Sterile



Zimmer® Biomet® Replacements

Compatible with Zimmer® Biomet® Knee Instrumentation Systems

Order Number	Replacement	Diameter (mm)	Length (mm)	Description	Packaging
ZB00-5901-020-00	00-5901-020-00	3.2	75.0	Headless Drill Pin (Hex End & Trocar Point)	4 Pins Per Pack Sterile

Biomet® is a registered trademark of Biomet, Inc., or its affiliates.

DePuy® and Synthes® are registered trademarks of DePuy Synthes, Inc., a Johnson & Johnson Company, or its affiliates.

Zimmer® is a registered trademark of Zimmer, Inc. or its affiliates.

Call 800.535.6638 to order today! Visit us online at BrasselerUSAMedical.com

