



**BRASSELER**  
USA® SURGICAL INSTRUMENTATION

# BRASSELER® NAVIGATION & EX-PINS

Brasseler USA Surgical Instrumentation is pleased to add new Navigation and Ex-Pins to our Stryker® Replacement Ortho Accessories line with ST6007-003-100, ST6007-003-150, ST6007-103-110, and ST6007-104-150. These replacement Navigation and Ex-Pins are compatible with Stryker® OrthoLock™ Navigation System.

Navigation pins are used to anchor the OrthoLock™ device on the femur and tibia bones for a TKA procedure and on the crest of the pelvis bone for a THA procedure. The patient trackers are then positioned on the OrthoLock™ device.

## ► COMPATIBLE WITH STRYKER® ORTHOLOCK™ NAVIGATION SYSTEM

Diameter	Length	Order Number	OEM Catalog Number
mm	mm	Sterile (2 Per Pack)	Stryker®
3	100	ST6007-003-100	6007-003-100
3	150	ST6007-003-150	6007-003-150
3	110	ST6007-103-110	6007-103-110
4	150	ST6007-104-150	6007-104-150

*Coming Soon!*

*Stryker® and OrthoLock™ are trademarks of Stryker Corporation or its affiliates.*

Call **800.535.6638** to order today!  
Visit us online at [BrasselerUSAMedical.com](http://BrasselerUSAMedical.com)



©2021 Brasseler USA. All rights reserved. Printed in the U.S.A.

BM-5186-01.21-2M-DXM

