



BRASELER
USA[®] SURGICAL INSTRUMENTATION

**MADE IN
USA** ★

NEW! **TELESCOPING BURS**

More ways to save! Introducing Brasseler Surgical Instrumentation's comprehensive line of high-speed telescoping burs.

- ▶ Compatible with Stryker[®] Signature Portfolio Elite Attachments
- ▶ Manufactured in the USA from high-quality tool steel and carbide
- ▶ Telescoping shaft with up to 5 adjustable exposure levels
- ▶ Featuring multiple head types including: Round Fluted, Neuro Match Head, and Cross-Cut Round

Make the switch today! Our seamless and easy conversion process makes reducing your bottom line simple. Feel confident making the switch without sacrificing quality or performance. **SAVE UP TO 60% vs. OEM!***

*Data on file.

- ▶ Visit our website at BrasselerUSAMedical.com
To order call **800.535.6638**



©2022 Brasseler U.S.A. Medical, LLC. All rights reserved. Printed in the U.S.A.



**Switch to
Brasseler &
Save!***

BM-5536-01.22

▶ BRASSLER USA TELESCOPING BURS

Compatible with Stryker® Signature Portfolio Elite Attachments



Brasseler® Part No.	Stryker® Part No.	Description	Head Diameter (mm)
BR5820-010-120	5820-010-120	Round Fluted Bur – Less Aggressive	2.0
BR5820-010-130	5820-010-130	Round Fluted Bur – Less Aggressive	3.0
BR5820-010-140	5820-010-140	Round Fluted Bur – Less Aggressive	4.0
BR5820-010-150	5820-010-150	Round Fluted Bur – Less Aggressive	5.0
BR5820-010-160	5820-010-160	Round Fluted Bur – Less Aggressive	6.0
BR5820-010-170	5820-010-170	Round Fluted Bur – Less Aggressive	7.0
BR5820-010-230	5820-010-230	Round Fluted Bur – Aggressive	3.0
BR5820-010-240	5820-010-240	Round Fluted Bur – Aggressive	4.0
BR5820-010-250	5820-010-250	Round Fluted Bur – Aggressive	5.0
BR5820-107-020	5820-107-020	Neuro Match Head Drill	2.0
BR5820-107-025	5820-107-025	Neuro Match Head Drill	2.5
BR5820-107-030	5820-107-030	Neuro Match Head Drill	3.0
BR5820-107-040	5820-107-040	Neuro Match Head Drill	4.0
BR5820-107-430	5820-107-430	Neuro Match Head Drill – Less Aggressive	3.0
BR5820-107-430S1	5820-107-430S1	Neuro Match Head Drill – Specialty	3.0
BR5820-107-530	5820-107-530	Neuro Match Head Drill – Exact	3.0
BR5820-107-530C	5820-107-530C	Neuro Match Head Drill – Carbide	3.0
BR5820-010-640	5820-010-640	Cross-Cut Round Bur	4.0
BR5820-010-650	5820-010-650	Cross-Cut Round Bur	5.0
BR5820-010-660	5820-010-660	Cross-Cut Round Bur	6.0

Stryker® is a registered trademark of Stryker Corporation or its affiliates.

▶ Visit our website at BrasselerUSAMedical.com
To order call **800.535.6638**

 **BRASSLER**
USA® SURGICAL INSTRUMENTATION