



**BRASSELER**  
USA<sup>®</sup> SURGICAL

## **NEW!** **TELESCOPING BURS**

More ways to save! Introducing Brasseler USA Surgical's comprehensive line of high-speed telescoping burs.

- ▶ Compatible with Stryker<sup>®</sup> Signature Portfolio Elite Attachments
- ▶ Manufactured in the USA from high-quality tool steel and carbide
- ▶ Telescoping shaft with up to 5 adjustable exposure levels
- ▶ Featuring multiple head types including: Round Fluted, Neuro Match Head, Cross-Cut and Diamond

Make the switch today! Our seamless and easy conversion process makes reducing your bottom line simple. Because Brasseler products provide exceptional performance and reliability, you can feel confident making the switch. **SAVE UP TO 60% vs. OEM!**\*

\*Based on analysis of publicly available retail pricing.  
Data on file and available upon request.


**MADE IN**  
**USA** 



**Switch to  
Brasseler &  
Save!\***

# Brasseler USA® Telescoping Burs

Compatible with Stryker® Signature Portfolio Elite Attachments



Brasseler® Part No.	Compare to Stryker® Part No.	Description	Head Diameter (mm)	Head Length (mm)
<b>BR5820-010-120</b>	5820-010-120	Round Fluted Bur – Less Aggressive	2.0	2.0
<b>BR5820-010-130</b>	5820-010-130	Round Fluted Bur – Less Aggressive	3.0	3.0
<b>BR5820-010-140</b>	5820-010-140	Round Fluted Bur – Less Aggressive	4.0	4.0
<b>BR5820-010-150</b>	5820-010-150	Round Fluted Bur – Less Aggressive	5.0	5.0
<b>BR5820-010-160</b>	5820-010-160	Round Fluted Bur – Less Aggressive	6.0	6.0
<b>BR5820-010-170</b>	5820-010-170	Round Fluted Bur – Less Aggressive	7.0	7.0
<b>BR5820-010-230</b>	5820-010-230	Round Fluted Bur – Aggressive	3.0	3.0
<b>BR5820-010-240</b>	5820-010-240	Round Fluted Bur – Aggressive	4.0	4.0
<b>BR5820-010-250</b>	5820-010-250	Round Fluted Bur – Aggressive	5.0	5.0
<b>BR5820-107-020</b>	5820-107-020	Neuro Match Head Drill	2.0	3.3
<b>BR5820-107-025</b>	5820-107-025	Neuro Match Head Drill	2.5	3.7
<b>BR5820-107-030</b>	5820-107-030	Neuro Match Head Drill	3.0	3.8
<b>BR5820-107-040</b>	5820-107-040	Neuro Match Head Drill	4.0	5.1
<b>BR5820-107-430</b>	5820-107-430	Neuro Match Head Drill – Less Aggressive	3.0	3.8
<b>BR5820-107-430S1</b>	5820-107-430S1	Neuro Match Head Drill – Specialty	3.0	3.8
<b>BR5820-107-530</b>	5820-107-530	Neuro Match Head Drill – Exact	3.0	3.8
<b>BR5820-107-530C</b>	5820-107-530C	Neuro Match Head Drill – Carbide	3.0	3.8
<b>BR5820-010-640</b>	5820-010-640	Cross-Cut Round Bur	4.0	4.0

Brasseler® Part No.	Compare to Stryker® Part No.	Description	Head Diameter (mm)	Head Length (mm)
<b>BR5820-010-650</b>	5820-010-650	Cross-Cut Round Bur	5.0	5.0
<b>BR5820-010-660</b>	5820-010-660	Cross-Cut Round Bur	6.0	6.0
<b>BR5820-012-020</b>	5820-012-020	Round Diamond Bur	2.0	2.0
<b>BR5820-012-030</b>	5820-012-030	Round Diamond Bur	3.0	3.0
<b>BR5820-012-040</b>	5820-012-040	Round Diamond Bur	4.0	4.0
<b>BR5820-012-050</b>	5820-012-050	Round Diamond Bur	5.0	5.0
<b>BR5820-012-060</b>	5820-012-060	Round Diamond Bur	6.0	6.0
<b>BR5820-013-020</b>	5820-013-020	Round Diamond Bur - Coarse	2.0	2.0
<b>BR5820-013-030</b>	5820-013-030	Round Diamond Bur - Coarse	3.0	3.0
<b>BR5820-013-040</b>	5820-013-040	Round Diamond Bur - Coarse	4.0	4.0
<b>BR5820-013-050</b>	5820-013-050	Round Diamond Bur - Coarse	5.0	5.0
<b>BR5820-013-140</b>	5820-013-140	Round Diamond Bur - X-Coarse	4.0	4.0
<b>BR5820-013-150</b>	5820-013-150	Round Diamond Bur - X-Coarse	5.0	5.0
<b>BR5820-012-310</b>	5820-012-310	Tapered Round Diamond Bur	1.0	1.0
<b>BR5820-012-315</b>	5820-012-315	Tapered Round Diamond Bur	1.5	1.5
<b>BR5820-012-320</b>	5820-012-320	Tapered Round Diamond Bur	2.0	2.0
<b>BR5820-107-130</b>	5820-107-130	Neuro Diamond Drill	3.0	3.8
<b>BR5820-066-216</b>	5820-066-216	Cross Cut Fissure Bur - Tungsten Carbide	1.6	6.2

Stryker® is a registered trademark of Stryker Corporation or its affiliates.

# MANAGE COSTS WITHOUT SACRIFICING QUALITY

- ✓ Achieve up to 40-60% Savings vs. OEM
- ✓ Compatible With Major OEM Power Systems
- ✓ Seamless & Easy Conversion Process
- ✓ Over 35 Years of Sales and Manufacturing Experience
- ✓ No Minimum Order Quantity
- ✓ Extensive Cross Referencing System by OEM
- ✓ Free 2nd Day Shipping on Orders Over \$100
- ✓ Free Evaluation Samples



Visit our website at [BrasselerUSAMedical.com](https://BrasselerUSAMedical.com)

To order call 800.535.6638 or fax 912.921.7578.

