



BRASSELER
USA[®] SURGICAL

MADE IN
USA ★

REPLACEMENT GUIDE WIRES

Brasseler USA Surgical is pleased to offer two partially threaded guide wires in the DePuy Synthes[®] compatible orthopedic accessories line. These guide wires are used to properly align screws and nails in orthopedic trauma procedures and are designed to fit in OEM procedural trays.

► **Compatible with DePuy Synthes[®] 3.5mm & 4.0mm Cannulated Screws**

- Fabricated from 17-4 stainless steel
- Diamond Point

Order Number	OEM Catalog Number	Diameter	Length
Sterile (1 Per Pack)	Depuy Synthes [®]	mm	mm
BR900-722	900.722	1.25	150.0

► **Compatible with DePuy Synthes[®] Proximal Femoral Nailing Systems**

- Fabricated from 316L Stainless Steel
- Self-Drilling Trocar Point

Order Number	OEM Catalog Number	Diameter	Length
Sterile (1 Per Pack)	Depuy Synthes [®]	mm	mm
BR357-399	357.399	3.2	400.0

DePuy Synthes[®] is a registered trademark of DePuy Synthes, Inc., a Johnson & Johnson Company, or its affiliates.

Call **800.535.6638** to order today!
Visit us online at BrasselerUSAMedical.com

