



BRASSELER
USA® SURGICAL

VersaDriver

PNEUMATIC OSTEOTOME SYSTEM



Clean, Fast Cement Removal for
the Most Challenging Surgeries

► FEATURES

- Knee & hip applications
- Cemented & non-cemented implants
- Precise cutting control
- Single or multi-strike cutting actions
- Blade limited to 3mm stroke
- Up to 3900 strikes/minute
- Hand piece & six Osteotome blades

MADE IN
USA

Distributed by Brasseler USA Surgical

Call 800.535.6638 to order today!

Visit us online at BrasselerSurgical.com






BRASSELER
USA® SURGICAL

VersaDriver

PNEUMATIC OSTEOTOME SYSTEM



PRODUCT CODE	ITEM DESCRIPTION
10-100	VersaDriver Pneumatic Osteotome System
10-001	VersaDriver Hand Piece
10-092	VersaDriver Case
10-093	VersaDriver Tray
10-095	Hose
10-069	Blade Extraction Slap Hammer
10-101	Flat Tip Blade 11mm x 50mm
10-102	Flat Tip Blade 15.5mm x 100mm
10-103	Round Tip Blade 11mm x 100mm
10-104	Round Tip Blade 15.5mm x 100mm
10-105	Flat Tip Blade 11mm x 100mm
10-106	Lateral Hip 18° Curved Tip Blade 11mm x 100mm
10-107	Medial Hip 18° Reverse Curved Tip Blade 11mm x 100mm
10-108	Offset Knee 20mm Round Tip Blade 5mm x 100mm
10-109	Gouge Curved Round Tip Blade 12mm Diam x 100mm
10-110	Gouge Curved Round Tip Blade 14mm Diam x 100mm
10-111	Cement Removal V Tip Blade 8mm x 100mm
10-112	Cement Removal Parallel V Tip Blade 15.5mm x 100mm