




BRASSELER
USA[®] SURGICAL

K-WIRES & STEINMANN PINS

Brasseler USA is proud to offer a comprehensive, high-quality K-Wire and Steinmann Pin product line. Fabricated from precision ground 316L implant grade stainless steel (ASTM F 138), our wires and pins are manufactured in the USA to very rigid specifications and subjected to strenuous inspection programs to ensure predictable, positive results. We are pleased to provide only the best quality products at an exceptional value.


- ▶ **Precision Manufactured for Easier Penetration**
- ▶ **Non-sterile (6/pack) & Sterile (1/pack)**
- ▶ **Multiple Point Styles & Sizes**

MADE IN
USA 




K-WIRES

SMOOTH K-WIRES							
Diameter		4" Length		6" Length		9" Length	
mm	in	Non-Sterile (6/pk)	Sterile (1/pk)	Non-Sterile (6/pk)	Sterile (1/pk)	Non-Sterile (6/pk)	Sterile (1/pk)
Single Trocar							
0.7	.028	KM172-14-28	KM172-14-28S	KM172-16-28	KM172-16-28S	N/A	N/A
0.9	.035	KM172-14-35	KM172-14-35S	KM172-16-35	KM172-16-35S	KM172-19-35	KM172-19-35S
1.1	.045	KM172-14-45	KM172-14-45S	KM172-16-45	KM172-16-45S	KM172-19-45	KM172-19-45S
1.4	.054	KM172-14-54	KM172-14-54S	KM172-16-54	KM172-16-54S	KM172-19-54	KM172-19-54S
1.6	.062	KM172-14-62	KM172-14-62S	KM172-16-62	KM172-16-62S	KM172-19-62	KM172-19-62S
Double Trocar							
0.7	.028	KM172-24-28	KM172-24-28S	KM172-26-28	KM172-26-28S	KM172-29-28	KM172-29-28S
0.9	.035	KM172-24-35	KM172-24-35S	KM172-26-35	KM172-26-35S	KM172-29-35	KM172-29-35S
1.1	.045	KM172-24-45	KM172-24-45S	KM172-26-45	KM172-26-45S	KM172-29-45	KM172-29-45S
1.4	.054	KM172-24-54	KM172-24-54S	KM172-26-54	KM172-26-54S	KM172-29-54	KM172-29-54S
1.6	.062	KM172-24-62	KM172-24-62S	KM172-26-62	KM172-26-62S	KM172-29-62	KM172-29-62S
Single Diamond							
0.7	.028	N/A	N/A	KM172-36-28	N/A	KM172-39-28	N/A
0.9	.035	N/A	N/A	KM172-36-35	N/A	KM172-39-35	N/A
1.1	.045	N/A	N/A	KM172-36-45	N/A	KM172-39-45	N/A
1.4	.054	N/A	N/A	N/A	N/A	KM172-39-54	N/A
1.6	.062	N/A	N/A	KM172-36-62	N/A	KM172-39-62	N/A
Double Diamond							
0.7	.028	KM172-44-28	KM172-44-28S	KM172-46-28	KM172-46-28S	KM172-49-28	KM172-49-28S
0.9	.035	KM172-44-35	KM172-44-35S	KM172-46-35	KM172-46-35S	KM172-49-35	KM172-49-35S
1.1	.045	KM172-44-45	KM172-44-45S	KM172-46-45	KM172-46-45S	KM172-49-45	KM172-49-45S
1.4	.054	KM172-44-54	KM172-44-54S	KM172-46-54	KM172-46-54S	KM172-49-54	KM172-49-54S
1.6	.062	KM172-44-62	KM172-44-62S	KM172-46-62	KM172-46-62S	KM172-49-62	KM172-49-62S
Extra Sharp, Double Diamond							
0.7	.028	N/A	KM172-44-28ESS	N/A	N/A	N/A	N/A
0.9	.035	N/A	KM172-44-35ESS	N/A	N/A	N/A	N/A
1.1	.045	N/A	KM172-44-45ESS	N/A	KM172-46-45ESS	N/A	N/A
1.6	.062	N/A	KM172-44-62ESS	N/A	KM172-46-62ESS	N/A	N/A




Single Trocar

3-Sided Point; Point on One End




Double Trocar

3-Sided Point; Point on Both Ends



Single Diamond

2-Sided Point; Point on One End



Double Diamond

2-Sided Point; Point on Both Ends



Double Diamond, Extra Sharp


2-Sided Point; Point on Both Ends

STEINMANN PINS

FULLY THREADED K-WIRES			
Diameter		9" Length	
mm	in	Non-Sterile (6/pk)	Sterile (1/pk)
Single Trocar			
0.9	.035	KM173-19-35	KM173-19-35S
1.1	.045	KM173-19-45	KM173-19-45S
1.4	.054	N/A	N/A
1.6	.062	KM173-19-62	KM173-19-62S
Double Trocar			
0.9	.035	KM173-29-35	KM173-29-35S
1.1	.045	KM173-29-45	KM173-29-45S
1.4	.054	KM173-29-54	KM173-29-54S
1.6	.062	KM173-29-62	KM173-29-62S


FULLY THREADED K-WIRES			
Diameter		9" Length	
mm	in	Non-Sterile (6/pk)	Sterile (1/pk)
Single Diamond			
0.9	.035	KM173-39-35	KM173-39-35S
1.1	.045	KM173-39-45	KM173-39-45S
1.4	.054	KM173-39-54	KM173-39-54S
1.6	.062	KM173-39-62	KM173-39-62S
Double Diamond			
0.9	.035	KM173-49-35	KM173-49-35S
1.1	.045	KM173-49-45	KM173-49-45S
1.4	.054	KM173-49-54	KM173-49-54S
1.6	.062	KM173-49-62	KM173-49-62S

SMOOTH STEINMANN PINS						
Diameter		9" Length		12" Length		
mm	in	Non-Sterile (6/pk)	Sterile (1/pk)	Sterile (1/pk)	OEM Part Number / Zimmer® Biomet®	
Single Trocar						
2.0	5/64"	KM168-19-56	KM168-19-56S	—	—	
2.4	3/32"	KM168-19-33	KM168-19-33S	—	—	
2.8	7/64"	KM168-19-76	KM168-19-76S	—	—	
3.2	1/8"	KM168-19-18	KM168-19-18S	—	—	
3.6	9/64"	KM168-19-96	KM168-19-96S	—	—	
4.0	5/32"	KM168-19-53	KM168-19-53S	—	—	
4.8	3/16"	KM168-19-31	KM168-19-31S	—	—	
Double Trocar						
2.0	5/64"	KM168-29-56	KM168-29-56S	—	—	
2.4	3/32"	KM168-29-33	KM168-29-33S	—	—	
2.8	7/64"	KM168-29-76	KM168-29-76S	—	—	
3.2	1/8"	KM168-29-18	KM168-29-18S	—	—	
3.6	9/64"	KM168-29-96	KM168-29-96S	—	—	
4.0	5/32"	KM168-29-53	KM168-29-53S	—	—	
4.8	3/16"	KM168-29-31	KM168-29-31S	—	—	
Single Diamond						
2.0	5/64"	KM168-39-56	KM168-39-56S	—	—	
2.4	3/32"	KM168-39-33	KM168-39-33S	—	—	
2.8	7/64"	KM168-39-76	KM168-39-76S	—	—	
3.2	1/8"	KM168-39-18	KM168-39-18S	—	—	
3.6	9/64"	KM168-39-96	KM168-39-96S	—	—	
4.0	5/32"	KM168-39-53	KM168-39-53S	—	—	
4.8	3/16"	KM168-39-31	KM168-39-31S	—	—	
Double Diamond						
2.0	5/64"	KM168-49-56	KM168-49-56S	—	—	
2.4	3/32"	KM168-49-33	KM168-49-33S	—	—	
2.8	7/64"	KM168-49-76	KM168-49-76S	—	—	
3.2	1/8"	KM168-49-18	KM168-49-18S	—	—	
3.6	9/64"	KM168-49-96	KM168-49-96S	—	—	
4.0	5/32"	KM168-49-53	KM168-49-53S	ZB27-362093	7-362093	
4.8	3/16"	KM168-49-31	KM168-49-31S	—	—	




Single Trocar

3-Sided Point; Point on One End




Double Trocar

3-Sided Point; Point on Both Ends



Single Diamond

2-Sided Point; Point on One End



Double Diamond

2-Sided Point; Point on Both Ends



MANAGE COSTS WITHOUT SACRIFICING QUALITY

- ✓ Achieve up to 40-60% Savings vs. OEM
- ✓ Compatible With Major OEM Power Systems
- ✓ Seamless & Easy Conversion Process
- ✓ Over 35 Years of Sales and Manufacturing Experience
- ✓ No Minimum Order Quantity
- ✓ Extensive Cross Referencing System by OEM
- ✓ Free Evaluation Samples



Visit our website at [BrasselerSurgical.com](https://www.BrasselerSurgical.com)

To order call 800.535.6638.

