



BRASSELER
USA® SURGICAL

QUICK CONNECT TWIST DRILLS

Brasseler USA Surgical, your orthopedic solution provider, is proud to offer an extensive line of replacement Quick Connect Twist Drills. Available in a solid or cannulated cross section, our Quick Connect Twist Drills are compatible with DePuy® Synthes® Quick Connect Couplings. We are pleased to provide only the best quality products at an exceptional value.

- ▶ Autoclavable Stainless Steel
- ▶ Sterile & Non-Sterile Options
- ▶ Significant Cost Savings
- ▶ One Piece Construction

Make the switch today and SAVE UP TO 60% vs. OEM!*

*Based on analysis of publicly available retail pricing.
Data on file and available upon request.

**MADE IN
USA** 



**Switch to
Brasseler &
Save!***



▶ Cannulated

Brasseler Part No.	OEM Part No.	Diameter (mm)	Cannulation (mm)	Total Length (mm)	Effective Length (mm)	Number of Flutes
Non-sterile (1/pk)	Non-Sterile					
KM166-310-65	310.65	3.2	1.7	170	139	3
KM166-310-66	310.66	4.5	1.7	170	139	3
KM166-310-67	310.67	2.7	1.4	160	130	3
KM166-310-68	310.68	3.5	1.4	160	130	3
KM166-310-221	310.221	2.0	1.2	150	50	3

▶ Solid

Brasseler Part No.		OEM Part No.		Diameter (mm)	Total Length (mm)	Effective Length (mm)	Number of Flutes
Non-sterile (1/pk)	Sterile (1/pk)	Non-Sterile	Sterile				
KM166-310-11	KM166-310-11S	310.11	310.11S	1.1	60	35	2
KM166-310-15	KM166-310-15S	310.15	310.15S	1.5	87	61	2
KM166-310-16	KM166-310-16S	310.16	310.16S	1.5	111	86	2
KM166-310-19	KM166-310-19S	310.19	310.19S	2	100	75	2
KM166-310-21	KM166-310-21S	310.21	310.21S	2	125	100	2
KM166-310-23	KM166-310-23S	310.23	310.23S	2.5	181	156	2
KM166-310-25	KM166-310-25S	310.25	310.25S	2.5	111	86	2
KM166-310-26	KM166-310-26S	310.26	310.26S	2.7	102	76	2
KM166-310-28	KM166-310-28S	310.28	310.28S	2.7	126	100	2
KM166-310-29	KM166-310-29S	310.29	310.29S	3.2	196	170	2
KM166-310-31	KM166-310-31S	310.31	310.31S	3.2	147	121	2
KM166-310-35	KM166-310-35S	310.35	310.35S	3.5	112	86	2
KM166-310-37	KM166-310-37S	310.37	310.37S	3.5	196	171	2
KM166-310-44	KM166-310-44S	310.44	310.44S	4.5	147	121	2
KM166-310-48	KM166-310-48S	310.48	310.48S	4.5	195	170	2
KM166-310-288	N/A	310.288	N/A	2.8	165	140	2
KM166-310-509	N/A	310.509	N/A	1.8	110	85	2
KM166-310-530	KM166-310-530S	310.53	310.530S	2.4	100	75	2
KM166-315-28	N/A	315.28	N/A	2.7	127	102	3
KM166-323-062	N/A	323.62	N/A	2	140	115	3
KM166-324-214	N/A	324.214	N/A	2.8	197	171	3
KM166-357-407	N/A	357.407	N/A	4	260	235	3

Call 800.535.6638 to order today!
 Visit us online at BrasselerSurgical.com

