



BRASSELER
USA® SURGICAL

THE SMART CHOICE

EZ X

ACETABULAR CUP REMOVAL SYSTEM

for use in hip revision surgery



▶ FEATURES

- Compatible with any powered handpiece
- Removes cemented and non-cemented cups
- 29 blade options: size 46-74 in 1mm increments



▶ BENEFITS

- Reduces surgical procedural time
- Minimizes bone loss
- Simplifies complex hip revision surgeries, improving patient recovery and mobility
- Precise cut with predictable results
- Reduces the risk of acetabular fractures





PRESERVATION & PROFICIENCY



▶ PRODUCT OFFERING

ITEM NAME	ORDER CODE	ITEM DESCRIPTION
Universal Shaft	BREZS101	Shaft to assemble with Side Handle and Blade
Side Handle	BREZH100	Handle to hold Universal Shaft and forward Blade
Caliper Tool	BRXST2000M	Tool to size the Blade and Liners
Blade	BRXTB460 ~740	Blade size 46mm - 74mm (Available in 1 mm increments)
Raised Liner	BRXTL37A ~ 65A	Raised Liner 37mm - 65mm (Available in 1 mm increments)
Finishing Liner	BRXTL37B ~ 65B	Finishing Liner 37mm - 65mm (Available in 1 mm increments)
Sterilization Tray	BRXZT1000	Sterilization Tray

▶ Visit our website at BrasselerSurgical.com
To order call 800.535.6638

