




BRASSELER
USA® SURGICAL

ORTHOPAEDIC PIN & DRILL PACKS

Fabricated from surgical grade stainless steel and subjected to strenuous inspection programs to assure exceptional results



Compatible with Stryker® Knee Instrumentation Systems

OEM Part Number	Brasseler Part Number	Diameter (mm)	Length (mm)	Description	Packaging
7650-2038 7650-1038	BMA-S90008S	3.1	88.9	Headless Pin (Square End/Trocar Point)	4 Pins Per Pack Sterile
7650-2038A 7650-1038A	BMA-S90010S	3.1	88.9	Headless Pin, Fluted (Square End/Trocar Point)	4 Pins Per Pack Sterile
6541-4-003 6541-004-003	BMA-S90012S	3.1	76.2	Headless Pin (Square End/Trocar Point)	4 Pins Per Pack Sterile

Compatible with Stryker® OrthoLock™ Navigation System

OEM Part Number	Brasseler Part Number	Diameter (mm)	Length (mm)	Description	Packaging
6007-003-100	ST6007-003-100	3	100	Navigation Pin (Square End/Trocar Point)	2 Pins Per Pack Sterile
6007-003-150	ST6007-003-150	3	150	Navigation Pin (Square End/Trocar Point)	2 Pins Per Pack Sterile
6007-103-110	ST6007-103-110	3	110	Navigation Pin (Square End/Trocar Point)	2 Pins Per Pack Sterile
6007-104-150	ST6007-104-150	4	150	Navigation Pin (Square End/Trocar Point)	2 Pins Per Pack Sterile

Stryker® and OrthoLock™ are trademarks of Stryker Corporation or its affiliates.

Call 800.535.6638 to order today! Visit us online at BrasselerUSAMedical.com





BRASSELER
USA® SURGICAL

ORTHOPAEDIC PIN & DRILL PACKS

Fabricated from surgical grade stainless steel and subjected to strenuous inspection programs to assure exceptional results

Compatible with DePuy® Synthes® Knee Instrumentation Systems

OEM Part Number	Brasseler Part Number	Diam. (mm)	Length (mm)	Description	Packaging
2490-95-000	BMA-S90006S	3.1	76.2	Fixation Pin (Round End/Trocar Point)	4 Pins Per Pack Sterile
86-4192-000	BMA-S90007S	3.1	76.2	Fixation Pin (Triangular End/Trocar Point) & Twist Drills (Round End)	4 Pins & 2 Drills Per Pack Sterile
86-9117	DP86-9117	3.1	76.2	Fixation Pin (Round End/Trocar Point)	10 Pins Per Pack Non-Sterile
2544-00-111	BR2544-00-111	3.2	67	Threaded Headed & Universal Pin	2 Threaded Headed & 4 Universal Pins Per Pack, Sterile
9505-02-302	BR9505-02-302	3.2	76.2	Threaded Pin	4 Pins Per Pack Sterile
9505-02-303	BR9505-02-303	3.2	67	Threaded Headed Pin	4 Pins Per Pack Sterile

Compatible with Zimmer® Biomet® Knee Instrumentation Systems

OEM Part Number	Brasseler Part Number	Diam. (mm)	Length (mm)	Description	Packaging
00-5901-020-00	ZB00-5901-020-00	3.2	75.0	Headless Drill Pin (Hex End/Trocar Point)	4 Pins Per Pack Sterile

Compatible with Smith+Nephew® GENESIS® Total Knee System

OEM Part Number	Brasseler Part Number	Diam. (mm)	Length (mm)	Description	Packaging
114968	BML-S90004	3.2	127/76	Twist Drill (Round End) & Fixation Pin (Round End/Cone Point)	1 Drill & 3 Pins Per Pack Sterile
71210002	BMA-S90006S	3.1	76.2	Fixation Pin (Round End/Trocar Point)	4 Pins Per Pack Sterile

Biomet® is a registered trademark of Biomet, Inc., or its affiliates.
DePuy® and Synthes® are registered trademarks of DePuy Synthes, Inc., a Johnson & Johnson Company, or its affiliates.
Smith+Nephew® and Genesis® are registered trademarks of Smith & Nephew, Inc., or its affiliates.
Zimmer® is a registered trademark of Zimmer, Inc. or its affiliates.

Call 800.535.6638 to order today! Visit us online at BrasselerUSAMedical.com

