



BRASSELER
USA® SURGICAL


**Switch to
Brasseler &
Save!***

SMALL BONE SAW BLADES

- ▶ *Manufactured in Savannah, GA*
- ▶ *Compatible with major OEM power systems*
- ▶ *Extensive cross referencing system by OEM*
- ▶ *Seamless and easy conversion process*

New Blades Now Available!

SAVE UP TO 60% vs. OEM



**MADE IN
USA**

*Based on analysis of publicly available retail pricing. Data on file upon request



Increasing Your Bottom Line By Reducing OR Costs

Our comprehensive line of high-quality replacement surgical small bone saw blades are compatible with most Original Equipment Manufacturers (OEMs)



Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1,
TPX®, & Command® 2

Compatible with Stryker® Reciprocating Saw Systems:
RemB®, CORE™, TPS™, F1, TPX®, & Command® 2





Compatible with Stryker® Reciprocating Saw Systems:
Micro

Compatible with ConMed® Linvatec® Sagittal Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; Brasseler USA®,
MicroAire®, & Smith+Nephew® Small Bone/Micro Sagittal Saw Systems



Compatible with ConMed® Linvatec® Oscillating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; Brasseler USA®,
MicroAire®, & Smith+Nephew® Small Bone/Micro Oscillating Saw Systems

Compatible with ConMed® Linvatec® Reciprocating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; Brasseler USA®,
MicroAire®, & Smith+Nephew® Small Bone/Micro Reciprocating Saw Systems



Compatible with ConMed® Linvatec® Mini-Driver™ K100 & K200 Sagittal Saw
Systems: K120 Attachment

Compatible with ConMed® Linvatec® Mini-Driver™ K100 & K200 Sagittal Saw
Systems: K220 Attachment



Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX®, & Command® 2



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)
2296-003-051	KM3-51	18.5	6.5	0.38	0.43
2296-003-101	KM3-101	10.0	9.0	0.38	0.43
2296-003-102	KM3-102	18.5	16.5	0.38	0.43
2296-003-103	KM3-103	15.0	7.0	0.38	0.43
2296-003-104	KM3-104	34.5	16.5	0.38	0.43
2296-003-105	KM3-105	18.5	9.0	0.38	0.43
2296-003-106	KM3-106	34.5	13.0	0.38	0.43
2296-003-108	KM3-108	39.0	13.0	0.38	0.43
2296-003-109S2	KM3-109S2	14.3	15.9	0.38	0.68
2296-003-111	KM3-111	25.0	9.0	0.38	0.43
2296-003-114	KM3-114	18.5	7.0	0.38	0.43
2296-003-115	KM3-115	29.5	7.0	0.38	0.43
2296-003-124	KM3-124	43.0	10.0	0.38	0.61
2296-003-125	KM3-125	31.0	9.0	0.38	0.43
2296-003-206	KM3-206	34.5	12.0	0.38	0.43
2296-003-212	KM3-212	18.0	5.5	0.38	0.43
2296-003-225	KM3-225	31.0	9.0	0.38	0.43
2296-003-254	KM3-254	41.0	19.5	0.38	0.68
2296-003-255	KM3-255	18.5	9.0	0.38	0.61
2296-003-410	KM3-410	11.5	5.5	0.38	0.43
2296-003-411	KM3-411	11.5	9.0	0.38	0.43
2296-003-412	KM3-412	18.0	5.5	0.38	0.43
2296-003-413	KM3-413	18.0	2.5	0.38	0.43
2296-003-414	KM3-414	25.0	5.5	0.38	0.43
2296-003-511	KM3-511	25.0	9.0	0.38	0.64
2296-003-515	KM3-515	42.0	13.4	0.38	0.61
2296-003-524	KM3-524	25.0	5.5	0.38	0.61
2296-003-525	KM3-525	31.0	9.0	0.38	0.64
2296-033-101	KM33-101	10.0	9.0	0.51	0.56
2296-033-102	KM33-102	18.5	16.5	0.51	0.56

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX®, & Command® 2



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)
2296-033-103	KM33-103	15.0	7.0	0.51	0.56
2296-033-104	KM33-104	34.5	16.5	0.51	0.56
2296-033-105	KM33-105	18.5	9.0	0.51	0.56
2296-033-106	KM33-106	34.5	13.0	0.51	0.56
2296-033-108	KM33-108	39.0	13.0	0.51	0.56
2296-033-111	KM33-111	25.0	9.0	0.51	0.56
2296-033-114	KM33-114	18.5	7.0	0.51	0.56
2296-033-115	KM33-115	29.5	7.0	0.51	0.56
2296-033-125	KM33-125	31.0	9.0	0.51	0.56
2296-033-206	KM33-206	34.5	12.0	0.51	0.56
2296-033-212	KM33-212	18.0	5.5	0.51	0.56
2296-033-225	KM33-225	31.0	9.0	0.51	0.56
2296-033-410	KM33-410	11.5	5.5	0.51	0.56
2296-033-411	KM33-411	11.5	9.0	0.51	0.56
2296-033-412	KM33-412	18.0	5.5	0.51	0.56
2296-033-413	KM33-413	18.0	2.5	0.51	0.56
2296-033-414	KM33-414	25.0	5.5	0.51	0.56
2296-033-520	KM33-520	35.0	9.0	0.63	0.86
2296-033-521	KM33-521	25.0	9.0	0.51	0.74
5400-003-105	BR5400-003-105	18.5	9.0	0.38	0.56
5400-003-412	BR5400-003-412	18.0	5.5	0.25	0.35
5400-003-111	BR5400-003-111	25.0	9.0	0.38	0.56
5400-003-414	BR5400-003-414	25.0	5.5	0.25	0.35
5400-003-125	BR5400-003-125	31.0	9.0	0.38	0.56
5400-033-527	BR5400-033-527	45.0	12.0	0.51	0.56
5400-134-279	BR5400-134-279	30.0	30.0	0.38	0.63
5400-134-282	BR5400-134-282	53.0	32.0	0.38	0.63
5400-134-415	BR5400-134-415	41.5	19.0	0.38	0.63

Compatible with Stryker® Reciprocating Saw Systems:
RemB®, CORE™, TPS™, F1, TPX®, & Command® 2



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)
5100-137-123	KM5100-137-123	27.0	6.4	0.38	0.51
5100-137-133	KM5100-137-133	21.8	6.4	0.38	0.51
5100-137-233	KM5100-137-233	27.4	6	0.38	0.51

OEM Part Number	Brasseler Part Number	Length (mm)	Max. Cut Width (mm)	Detail
5100-037-113	KM5100-37-113	11	5.5	Small Tear/Cross Cut Rasp
5100-037-114	KM5100-37-114	13	6.5	Large Tear/Cross Cut Rasp
5100-037-115	KM5100-37-115	11	5.5	Small Tear Rasp
5100-037-116	KM5100-37-116	13	6.5	Large Tear Rasp

Compatible with Stryker® Reciprocating Saw Systems:
Micro



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)
1675-133	KM-062	27.2	6	0.64	1.00

OEM Part Number	Brasseler Part Number	Length (mm)	Max. Cut Width (mm)	Detail
1675-116	KM-800	13	6.5	Large Tear Rasp
1675-115	KM-801	11	5.5	Small Tear Rasp
1675-113	KM-802	11	5.5	Small Tear/Cross Cut Rasp

Compatible with ConMed® Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; Brasseler USA®, MicroAire®, & Smith+Nephew® Small Bone/Micro Sagittal Saw Systems



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)	Detail
5023-133	KM-3100C	41.5	14.0	0.40	0.61	Straight Edge, Coarse
5023-138	KM-3101	25.7	9.5	0.40	0.60	Straight Edge, Fine
5023-143	KM-3101C	25.7	9.5	0.38	0.61	Straight Edge, Coarse
5023-139	KM-3102	25.7	4.5	0.40	0.60	Straight Edge, Fine
5044-145	KM-3104	16.1	9.5	0.40	0.63	Straight Edge, Fine
5023-144	KM-3105	16.1	4.5	0.40	0.63	Straight Edge, Fine
5023-132	KM-3106	25.7	5.5	0.40	0.61	Straight Edge, Coarse
5023-141	KM-3107	34.3	6.0	0.40	0.60	Straight Edge, Fine
5023-183	KM-3183	25.7	9.5	0.40	0.60	Straight Edge, Fine
5023-131	KM-3231	41.5	19.5	0.40	0.61	Straight Edge, Coarse
5023-137	KM-3237	25.7	14.0	0.40	0.61	Straight Edge, Coarse
5023-140	KM-3240	41.5	9.5	0.40	0.60	Straight Edge, Fine/Tapered

Compatible with ConMed® Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; Brasseler USA®, MicroAire®, & Smith+Nephew® Small Bone/Micro Oscillating Saw Systems



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)	Detail
5023-162	KM-3000	18.5	5.5	0.40	0.60	Straight Edge, Fine
5023-164	KM-3001	28.5	5.5	0.40	0.60	Straight Edge, Fine
5023-169	KM-3002	34.7	10.0	0.40	0.60	Straight Edge, Fine
N/A	KM-3003	22.0	5.5	0.40	0.60	Straight Edge, Fine
5023-165	KM-3004	27.9	9.5	0.40	0.60	Straight Edge, Fine
5023-163	KM-3005	18.5	10.0	0.40	0.60	Straight Edge, Fine
5023-160	KM-3006	12.0	5.5	0.40	0.60	Straight Edge, Fine
5023-166	KM-3008	28.5	18.5	0.40	0.60	Radial Edge, Fine
5023-161	KM-3261	12.0	10.0	0.40	0.60	Straight Edge, Fine
5023-165	KM-3265	28.5	10.0	0.40	0.60	Radial Edge, Fine
5023-170	KM-3170	12.0	9.0	0.40	0.60	Crescentic, Fine
5023-171	KM-3171	25.5	13.8	0.40	0.60	Crescentic, Fine
5023-172	KM-3172	32.0	18.0	0.40	0.60	Crescentic, Fine
5023-176	KM-3176	32.0	13.8	0.40	0.60	Crescentic, Fine

Compatible with ConMed® Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; Brasseler USA®, MicroAire®, & Smith+Nephew® Small Bone/Micro Reciprocating Saw Systems



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)	Detail
5053-220	KM-060	25.5	7.6	0.40	0.64	Flat Shank
5053-221	KM-070	38.0	7.6	0.40	0.64	Flat Shank
5053-254	KM-54	25.0	7.4	0.40	0.60	Round Shank

OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Detail
5053-080	KM-80	11.0	5.5	Small Tear Rasp - Round Shank
5053-081	KM-81	13.0	6.9	Large Tear Rasp - Round Shank
5053-082	KM-82	11.0	5.5	Small Tear/ Cross Cut Rasp - Round Shank
5053-083	KM-83	13.0	6.9	Large Tear/ Cross Cut Rasp - Round Shank

Compatible with ConMed® Linvatec® Mini-Driver™ K100 & K200 Sagittal Saw Systems - K120 Attachment



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)	Detail
K121	KM-21K	15.1	6.4	0.71	0.89	Fine
K122	KM-22K	31.0	6.4	0.71	0.89	Fine
K123	KM-23K	15.1	10.4	0.71	0.89	Fine
K124	KM-24K	31.0	10.4	0.71	0.89	Fine
K131	KM-31K	15.1	6.4	0.38	0.56	Fine
K133	KM-33K	15.1	10.4	0.38	0.56	Fine
K134	KM-34K	31.0	10.4	0.38	0.56	Fine

Compatible with ConMed® Linvatec® Mini-Driver™ K100 & K200 Sagittal Saw Systems - K220 Attachment



OEM Part Number	Brasseler Part Number	Length (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)	Detail
K141	KM-41K	14.2	6.4	0.71	0.97	Fine
K143	KM-43K	14.2	10.2	0.71	0.97	Fine
K144	KM-44K	30.0	10.2	0.71	0.97	Fine
K148	KM-48K	27.4	9.5	0.38	0.65	Fine
K151	KM-51K	14.2	6.4	0.38	0.65	Fine

ConMed®, Linvatec®, Micro100™, Micro E®, MicroChoice®, MicroPower®, and Mini-Driver™ are trademarks of CONMED Corporation or its affiliates. MicroAire® is a registered trademark of MicroAire Surgical Instruments, LLC, or its affiliates. Smith+Nephew® is a registered trademark of Smith & Nephew, Inc., or its affiliates.



MANAGE COSTS WITHOUT SACRIFICING QUALITY

- ✓ Achieve up to 40-60% Savings vs. OEM
- ✓ Compatible With Major OEM Power Systems
- ✓ Seamless & Easy Conversion Process
- ✓ Over 35 Years of Sales and Manufacturing Experience
- ✓ No Minimum Order Quantity
- ✓ Extensive Cross Referencing System by OEM
- ✓ Free 2nd Day Shipping on Orders Over \$100
- ✓ Free Evaluation Samples



Visit our website at [BrasselerUSA.com](https://www.BrasselerUSA.com)

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

