



BRASSELER
USA[®] SURGICAL

**Switch to
Brasseler &
Save!***

TELESCOPING BURS

More ways to save! Introducing Brasseler USA Surgical's comprehensive line of high-speed telescoping burs.


- ▶ Compatible with Stryker[®] Signature Portfolio Elite Attachments
- ▶ Manufactured in the USA from high-quality tool steel and carbide
- ▶ Telescoping shaft with up to 5 adjustable exposure levels
- ▶ Featuring multiple head types including: Round Fluted, Neuro Match Head, Cross-Cut and Diamond

Our seamless and easy conversion process makes reducing your bottom line simple. Because Brasseler products provide exceptional performance and reliability, you can feel confident making the switch.

Make the switch today and SAVE UP TO 60% vs. OEM!*


*Based on analysis of publicly available retail pricing.
Data on file and available upon request.

**MADE IN
USA** 




Brasseler USA® Telescoping Burs

Compatible with Stryker® Signature Portfolio Elite Attachments



Compare to Stryker® Part No.	Brasseler® Part No.	Description	Head Diameter (mm)	Head Length (mm)
5820-010-120	BR5820-010-120	Round Fluted Bur - Less Aggressive	2.0	2.0
5820-010-130	BR5820-010-130	Round Fluted Bur - Less Aggressive	3.0	3.0
5820-010-140	BR5820-010-140	Round Fluted Bur - Less Aggressive	4.0	4.0
5820-010-150	BR5820-010-150	Round Fluted Bur - Less Aggressive	5.0	5.0
5820-010-160	BR5820-010-160	Round Fluted Bur - Less Aggressive	6.0	6.0
5820-010-170	BR5820-010-170	Round Fluted Bur - Less Aggressive	7.0	7.0
5820-010-230	BR5820-010-230	Round Fluted Bur - Aggressive	3.0	3.0
5820-010-240	BR5820-010-240	Round Fluted Bur - Aggressive	4.0	4.0
5820-010-250	BR5820-010-250	Round Fluted Bur - Aggressive	5.0	5.0
5820-107-020	BR5820-107-020	Neuro Match Head Drill	2.0	3.3
5820-107-025	BR5820-107-025	Neuro Match Head Drill	2.5	3.7
5820-107-030	BR5820-107-030	Neuro Match Head Drill	3.0	3.8
5820-107-040	BR5820-107-040	Neuro Match Head Drill	4.0	5.1
5820-107-430	BR5820-107-430	Neuro Match Head Drill - Less Aggressive	3.0	3.8
5820-107-430S1	BR5820-107-430S1	Neuro Match Head Drill - Specialty	3.0	3.8
5820-107-530	BR5820-107-530	Neuro Match Head Drill - Exact	3.0	3.8
5820-107-530C	BR5820-107-530C	Neuro Match Head Drill - Carbide	3.0	3.8
5820-010-640	BR5820-010-640	Cross Cut Round Bur	4.0	4.0
5820-010-650	BR5820-010-650	Cross Cut Round Bur	5.0	5.0
5820-010-660	BR5820-010-660	Cross Cut Round Bur	6.0	6.0
5820-012-020	BR5820-012-020	Round Diamond Bur	2.0	2.0
5820-012-030	BR5820-012-030	Round Diamond Bur	3.0	3.0
5820-012-040	BR5820-012-040	Round Diamond Bur	4.0	4.0
5820-012-050	BR5820-012-050	Round Diamond Bur	5.0	5.0
5820-012-060	BR5820-012-060	Round Diamond Bur	6.0	6.0
5820-013-030	BR5820-013-030	Round Diamond Bur - Coarse	3.0	3.0
5820-013-040	BR5820-013-040	Round Diamond Bur - Coarse	4.0	4.0



Compare to Stryker® Part No.	Brasseler® Part No.	Description	Head Diameter (mm)	Head Length (mm)
5820-013-050	BR5820-013-050	Round Diamond Bur - Coarse	5.0	5.0
5820-013-140	BR5820-013-140	Round Diamond Bur - X-Coarse	4.0	4.0
5820-013-150	BR5820-013-150	Round Diamond Bur - X-Coarse	5.0	5.0
5820-012-310	BR5820-012-310	Tapered Round Diamond Bur	1.0	1.0
5820-012-315	BR5820-012-315	Tapered Round Diamond Bur	1.5	1.5
5820-012-320	BR5820-012-320	Tapered Round Diamond Bur	2.0	2.0
5820-107-130	BR5820-107-130	Neuro Diamond Drill	3.0	3.8
5820-066-216	BR5820-066-216	Cross Cut Fissure Bur - Tungsten Carbide	1.6	6.2
5820-080-021	BR5820-080-021	Helicoidal Rasp	2.1	19.6
5820-066-012	BR5820-066-012	Cross Cut Fissure Bur	1.2	4.4
5820-066-016	BR5820-066-016	Cross Cut Fissure Bur	1.6	4.7
5820-066-021	BR5820-066-021	Cross Cut Fissure Bur	2.1	2.1
5820-015-040	BR5820-015-040	Egg Bur	4.0	9.9
5820-015-050	BR5820-015-050	Egg Bur	5.0	9.9
5820-070-011	BR5820-070-011	Straight Router	1.1	6.0
5820-103-010	BR5820-103-010	Wire Pass Drill	1.0	6.7
5820-103-015	BR5820-103-015	Wire Pass Drill	1.5	6.7
5820-103-020	BR5820-103-020	Wire Pass Drill	2.0	6.7
5820-012-020D	BR5820-012-020D	Round Diamond Bur - Medium Grit	2.0	2.0
5820-012-030D	BR5820-012-030D	Round Diamond Bur - Medium Grit	3.0	3.0
5820-012-040D	BR5820-012-040D	Round Diamond Bur - Medium Grit	4.0	4.0
5820-012-050D	BR5820-012-050D	Round Diamond Bur - Medium Grit	5.0	5.0
5820-110-020C	BR5820-110-020C	Round Fluted Bur Carbide	2.0	2.0
5820-110-030C	BR5820-110-030C	Round Fluted Bur Carbide	3.0	3.0
5820-110-040C	BR5820-110-040C	Round Fluted Bur Carbide	4.0	4.0



MANAGE COSTS WITHOUT SACRIFICING QUALITY

- ✓ Achieve up to 40-60% Savings vs. OEM
- ✓ Compatible With Major OEM Power Systems
- ✓ Seamless & Easy Conversion Process
- ✓ Over 35 Years of Sales and Manufacturing Experience
- ✓ No Minimum Order Quantity
- ✓ Extensive Cross Referencing System by OEM
- ✓ Free Evaluation Samples



Visit us at [BrasselerSurgical.com](https://www.BrasselerSurgical.com)

To order call 800.535.6638.

