



**BRASELER**  
USA® SURGICAL

**Switch to  
Brasseler &  
Save!\***

## TELESCOPING BURS

More ways to save! Introducing Brasseler USA Surgical's comprehensive line of high-speed telescoping burs.


- ▶ Compatible with Stryker® Signature Portfolio Elite Attachments
- ▶ Manufactured in the USA from high-quality tool steel and carbide
- ▶ Telescoping shaft with up to 5 adjustable exposure levels
- ▶ Featuring multiple head types including: Round Fluted, Neuro Match Head, Cross-Cut and Diamond

Our seamless and easy conversion process makes reducing your bottom line simple. Because Brasseler products provide exceptional performance and reliability, you can feel confident making the switch.

***Make the switch today and SAVE UP TO 60% vs. OEM!\****


\*Based on analysis of publicly available retail pricing.  
Data on file and available upon request.

**MADE IN  
USA** ★



► TELESCOPING BURS

Compatible with Stryker® Signature Portfolio Elite Attachments



| Compare to Stryker® Part No. | Brasseler® Part No. | Description                              | Head Diameter (mm) | Head Length (mm) |
|------------------------------|---------------------|------------------------------------------|--------------------|------------------|
| 5820-010-120                 | BR5820-010-120      | Round Fluted Bur - Less Aggressive       | 2.0                | 2.0              |
| 5820-010-130                 | BR5820-010-130      | Round Fluted Bur - Less Aggressive       | 3.0                | 3.0              |
| 5820-010-140                 | BR5820-010-140      | Round Fluted Bur - Less Aggressive       | 4.0                | 4.0              |
| 5820-010-150                 | BR5820-010-150      | Round Fluted Bur - Less Aggressive       | 5.0                | 5.0              |
| 5820-010-160                 | BR5820-010-160      | Round Fluted Bur - Less Aggressive       | 6.0                | 6.0              |
| 5820-010-170                 | BR5820-010-170      | Round Fluted Bur - Less Aggressive       | 7.0                | 7.0              |
| 5820-010-230                 | BR5820-010-230      | Round Fluted Bur - Aggressive            | 3.0                | 3.0              |
| 5820-010-240                 | BR5820-010-240      | Round Fluted Bur - Aggressive            | 4.0                | 4.0              |
| 5820-010-250                 | BR5820-010-250      | Round Fluted Bur - Aggressive            | 5.0                | 5.0              |
| 5820-107-020                 | BR5820-107-020      | Neuro Match Head Drill                   | 2.0                | 3.3              |
| 5820-107-025                 | BR5820-107-025      | Neuro Match Head Drill                   | 2.5                | 3.7              |
| 5820-107-030                 | BR5820-107-030      | Neuro Match Head Drill                   | 3.0                | 3.8              |
| 5820-107-040                 | BR5820-107-040      | Neuro Match Head Drill                   | 4.0                | 5.1              |
| 5820-107-430                 | BR5820-107-430      | Neuro Match Head Drill - Less Aggressive | 3.0                | 3.8              |
| 5820-107-430S1               | BR5820-107-430S1    | Neuro Match Head Drill - Specialty       | 3.0                | 3.8              |
| 5820-107-530                 | BR5820-107-530      | Neuro Match Head Drill - Exact           | 3.0                | 3.8              |
| 5820-107-530C                | BR5820-107-530C     | Neuro Match Head Drill - Carbide         | 3.0                | 3.8              |
| 5820-010-640                 | BR5820-010-640      | Cross Cut Round Bur                      | 4.0                | 4.0              |
| 5820-010-650                 | BR5820-010-650      | Cross Cut Round Bur                      | 5.0                | 5.0              |
| 5820-010-660                 | BR5820-010-660      | Cross Cut Round Bur                      | 6.0                | 6.0              |
| 5820-012-020                 | BR5820-012-020      | Round Diamond Bur                        | 2.0                | 2.0              |
| 5820-012-030                 | BR5820-012-030      | Round Diamond Bur                        | 3.0                | 3.0              |
| 5820-012-040                 | BR5820-012-040      | Round Diamond Bur                        | 4.0                | 4.0              |
| 5820-012-050                 | BR5820-012-050      | Round Diamond Bur                        | 5.0                | 5.0              |
| 5820-012-060                 | BR5820-012-060      | Round Diamond Bur                        | 6.0                | 6.0              |
| 5820-013-020                 | BR5820-013-020      | Round Diamond Bur - Coarse               | 2.0                | 2.0              |
| 5820-013-030                 | BR5820-013-030      | Round Diamond Bur - Coarse               | 3.0                | 3.0              |
| 5820-013-040                 | BR5820-013-040      | Round Diamond Bur - Coarse               | 4.0                | 4.0              |
| 5820-013-050                 | BR5820-013-050      | Round Diamond Bur - Coarse               | 5.0                | 5.0              |

| Compare to Stryker® Part No. | Brasseler® Part No. | Description                      | Head Diameter (mm) | Head Length (mm) |
|------------------------------|---------------------|----------------------------------|--------------------|------------------|
| 5820-013-140                 | BR5820-013-140      | Round Diamond Bur - X-Coarse     | 4.0                | 4.0              |
| 5820-013-150                 | BR5820-013-150      | Round Diamond Bur - X-Coarse     | 5.0                | 5.0              |
| 5820-012-310                 | BR5820-012-310      | Tapered Round Diamond Bur        | 1.0                | 1.0              |
| 5820-012-315                 | BR5820-012-315      | Tapered Round Diamond Bur        | 1.5                | 1.5              |
| 5820-012-320                 | BR5820-012-320      | Tapered Round Diamond Bur        | 2.0                | 2.0              |
| 5820-107-130                 | BR5820-107-130      | Neuro Diamond Drill              | 3.0                | 3.8              |
| 5820-066-216                 | BR5820-066-216      | Cross Cut Fissure Bur - Extended | 1.6                | 6.2              |
| 5820-080-021                 | BR5820-080-021      | Helicoidal Rasp                  | 2.1                | 19.6             |
| 5820-066-012                 | BR5820-066-012      | Cross Cut Fissure Bur            | 1.2                | 4.4              |
| 5820-066-016                 | BR5820-066-016      | Cross Cut Fissure Bur            | 1.6                | 4.7              |
| 5820-066-021                 | BR5820-066-021      | Cross Cut Fissure Bur            | 2.1                | 2.1              |
| 5820-015-040                 | BR5820-015-040      | Egg Bur                          | 4.0                | 9.9              |
| 5820-015-050                 | BR5820-015-050      | Egg Bur                          | 5.0                | 9.9              |
| 5820-070-011                 | BR5820-070-011      | Straight Router                  | 1.1                | 6.0              |
| 5820-103-010                 | BR5820-103-010      | Wire Pass Drill                  | 1.0                | 6.7              |
| 5820-103-015                 | BR5820-103-015      | Wire Pass Drill                  | 1.5                | 6.7              |
| 5820-103-020                 | BR5820-103-020      | Wire Pass Drill                  | 2.0                | 6.7              |
| 5820-012-020D                | BR5820-012-020D     | Round Diamond Bur - Medium Grit  | 2.0                | 2.0              |
| 5820-012-030D                | BR5820-012-030D     | Round Diamond Bur - Medium Grit  | 3.0                | 3.0              |
| 5820-012-040D                | BR5820-012-040D     | Round Diamond Bur - Medium Grit  | 4.0                | 4.0              |
| 5820-012-050D                | BR5820-012-050D     | Round Diamond Bur - Medium Grit  | 5.0                | 5.0              |
| 5820-110-020C                | BR5820-110-020C     | Round Fluted Bur Carbide         | 2.0                | 2.0              |
| 5820-110-030C                | BR5820-110-030C     | Round Fluted Bur Carbide         | 3.0                | 3.0              |
| 5820-110-040C                | BR5820-110-040C     | Round Fluted Bur Carbide         | 4.0                | 4.0              |
| 5820-110-050C                | BR5820-110-050C     | Round Fluted Bur Carbide         | 5.0                | 5.0              |
| 5820-110-060C                | BR5820-110-060C     | Round Fluted Bur Carbide         | 6.0                | 6.0              |



# MANAGE COSTS WITHOUT SACRIFICING QUALITY

- ✓ Achieve up to 40-60% Savings vs. OEM
- ✓ Compatible With Major OEM Power Systems
- ✓ Seamless & Easy Conversion Process
- ✓ Over 35 Years of Sales and Manufacturing Experience
- ✓ No Minimum Order Quantity
- ✓ Extensive Cross Referencing System by OEM
- ✓ Free Evaluation Samples



Visit us at [BrasselerSurgical.com](https://BrasselerSurgical.com)

To order call 800.535.6638.

