





BRASSELER
USA[®] SURGICAL



STERILIZATION TRAYS THAT KEEP YOU READY

Brasseler Surgical's
Sterilization Trays keep
instruments secure,
organized, and ready for
use in the OR.



Call **800.535.6638** to order today!
Visit us online at [BrasselerSurgical.com](https://www.BrasselerSurgical.com)

©2025 Brasseler USA. All rights reserved. Printed in the U.S.A.

NEW!


BRASSELER
 USA[®] SURGICAL

Quick Connect Drill Tray

Included in the Loaded Tray

Brasseler Product No.	Length (mm)	OEM Product No.	Quantity (ea)	Diameter (mm)
KM166-310-11	60	310.11	2	1.1
KM166-310-15	87	310.15	2	1.5
KM166-310-509	110	310.509	2	1.8
KM166-310-19	100	310.19	4	2.0
KM166-310-530	100	310.530	2	2.4
KM166-310-25	111	310.25	4	2.5
KM166-310-26	102	310.26	2	2.7
KM166-310-288	165	310.288	4	2.8
KM166-310-31	147	310.31	4	3.2
KM166-310-35	112	310.35	2	3.5
KM166-310-431	180	310.431	2	4.3
KM166-310-44T	152	310.44	2	4.5

Product No.	Item
BR-9708	QC Tray (Empty)
BR-9708K	QC Tray (Loaded)

- ▶ Holds 32 Quick Connect drills bits in 12 different sizes
- ▶ Available empty or loaded
- ▶ Drill bits may be interchanged with other available diameters, up to 197mm






Product No.	Item
BR-9704	Twist Drill Tray (Empty)
BR-9704K	Twist Drill Tray (Loaded)

TWIST DRILL TRAY

Included in the Loaded Tray


Product No. (Non-Sterile)	Product No. (Sterile)	Length (in)	Quantity (ea)	Diameter (mm)	Diameter (in)
KM166-00-11	KM166-10-11	4	2	1.1	3/64
KM166-00-00	KM166-10-00	5	2	1.6	1/16
KM166-00-01	KM166-10-01	5	4	2.0	5/64
KM166-00-02	KM166-10-02	5	4	2.4	3/32
KM166-00-03	KM166-10-03	5	4	2.8	7/64
KM166-00-04	KM166-10-04	5	2	3.2	1/8
KM166-00-05	KM166-10-05	5	2	3.6	9/64
KM166-00-06	KM166-10-06	5	2	4.0	5/32
KM166-00-45	KM166-10-45	5	2	4.5	11/64
KM166-00-07	KM166-10-07	5	2	4.8	3/16
KM166-00-60	KM166-10-60	5	1	6.0	15/64

- ▶ Holds 27 Twist Drills in 11 different sizes
- ▶ Available empty or loaded
- ▶ Compact design



K-Wire Sizing Options			
Length (in)	Quantity (ea)	Diameter (mm)	Diameter (in)
4	6	0.7	0.028
4	6	0.9	0.035
4	6	1.1	0.045
4	6	1.4	0.054
4	6	1.6	0.062
6	6	0.7	0.028
6	6	0.9	0.035
6	6	1.1	0.045
6	6	1.4	0.054
6	6	1.6	0.062

- ▶ Holds 60 K-Wires in 4 & 6 inch lengths
- ▶ Multiple product configurations
- ▶ Lift tab design for easy product access
- ▶ Bi-fold design with double-stacked holder





K-WIRE & STEINMANN PIN TRAYS

Product No.	Item
BR-9705	K-Wire Tray
BR-9706	Steinmann Pin Tray

*All trays sold unloaded


- ▶ Holds 90 K-Wires in 4, 6 & 9 inch lengths
- ▶ Holds 42 Steinmann Pins in 9 inch length
- ▶ Lift tab design for easy product access
- ▶ Bi-fold design with triple-stacked holder

K-Wire Tray Sizing Options

Length (in)	Quantity (ea)	Diameter (mm)	Diameter (in)
4, 6, 9	6	0.7	0.028
	6	0.9	0.035
	6	1.1	0.045
	6	1.4	0.054
	6	1.6	0.062

Steinmann Pin Tray Sizing Options

Length (in)	Quantity (ea)	Diameter (mm)	Diameter (in)
9	6	2.0	5/64
	6	2.4	3/32
	6	2.8	7/64
	6	3.2	1/8
	6	3.6	9/64
	6	4.0	5/32
	6	4.8	3/16





VALUE IN THE O.R. SAVINGS ON THE BOTTOM LINE.

- ✓ Achieve up to 40-60% Savings vs. OEM*
- ✓ Compatible With Major OEM Power Systems
- ✓ Seamless & Easy Conversion Process
- ✓ Over 35 Years of Sales & Manufacturing Experience
- ✓ No Minimum Order Quantity
- ✓ Extensive Cross Reference System by OEM

Call **800.535.6638** to order today!
Visit us online at BrasselerSurgical.com

*Based on analysis of publicly available retail pricing on. Data on file and available upon request.
©2025 Brasseler USA. All rights reserved. Printed in the U.S.A.



BM-5946 2025-09-15